Otto Bock®

QUALITY FOR LIFE

6Y40, 6Y41, 6Y42, 6Y70, 6Y71, 6Y75, 6Y80, 6Y81, 6Y85, 88L, 6Y90, 6Y92, 6Y93, 6Y60, 6Y61, 6Y510, 6Y511, 6Y520, 6Y521, 6Y522, 6Y523, 6Y430, 6Y400, 6Y416, 6Y414, 6Y540

Ⓓ	Otto Bock Liner	2
ⒼⒷ	Otto Bock Liners	7
Ⓕ	Liner Otto Bock	12
Ⓘ	Liner Otto Bock	18
Ⓔ	Liners de Otto Bock	24
ⓃⓁ	Otto Bock Liners	29
Ⓐ	Otto Bock Liner	35
Ⓝ	Otto Bock Liner	40
Ⓟ	Leje firmy Otto Bock	46
Ⓗ	Otto Bock Linerek	51
ⒸⓏ	Linery Otto Bock	57
Ⓡⓞ	Linere Otto Bock	63
ⓉⓇ	Otto Bock Liner	69
ⒼⓇ	Θήκες Otto Bock	74
Ⓡⓤⓢ	Культеприемники пр-ва Otto Bock	80
Ⓜ	奥托博克舒适型内衬套	86

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1 Intended Use

English

All Otto Bock liners are to be used for prosthetic fittings of transtibial amputation, except the TF Adapt Liner 6Y80/6Y81/6Y85 and the Akquire 6Y540. The 6Y80/6Y81/6Y85 and 6Y540 liners are to be used for prosthetic fittings of transfemoral amputations. Custom made liners are to be used according to their indication.

Note:

Liners that are provided with the Sanitized® silver antibacterial active substance show discolouration in 1~4% of patients. This discolouration can occur when the sweat contains sulphur-containing substances. The sulphur may appear in the sweat due to the consumption of the following foods:

- sulphur containing foods (such as protein products, cabbage, nuts or legume)
- processed foods preserved by sulfites (e.g. some dried fruits)

The discolouration occurs due to reaction of silver in the Sanitized® active substance with sulphur in the above foods. However, this is a completely natural reaction and can also be observed on jewellery made of silver or silver cutlery.

The discolouration has no influence of the liner, however, it cannot be removed by washing.

No harmful effects on residual limbs are known.

Due to this fact no returns based on discolouration are possible.

2 General Safety Instructions

Warning!

Keep the product away from children, since it presents a danger of suffocation!

Make sure not to wear the liner on injured skin!

Notice:

Attention! Otto Bock liners perform best when first use is before the printed "best before date".

3 Description

3.1 Silicone

Silicone liners 6Y40/41/42, SiliconeGel liners 6Y70/71/75, TF Adapt liner 6Y80/6Y81/6Y85, Custom made silicone liner: these liners are made of special silicones, which offer shock absorption and exceptional durability. Silicone non-locking liners can be used as compressive aid to reduce post-surgical oedema. Refer to your prescribing physician for treatment plan.

3.2 Copolymer

Copolymer liners 6Y90/92/93: these liners are made of a special TPE-Gel offering high reduction of shear forces and medical grade oil to create a skin-friendly and skin moisturising environment.

3.3 Polyurethane

Polyurethane (PUR) liners Profile 6Y510/511, Simplicity 6Y520/521/522/523, TechnoGel 6Y60/61, Akquire 6Y540 and Custom Liners 6Y400/414/416: these liners are made from special PUR Gel which offers exceptional shock absorption and yielding (floating) effect to provide a protective environment for limbs.

Note:

Liners normally feel softer after use than the newly received ones. This applies to all three liner materials.

4 Indication

4.1 Silicone and SiliconeGel

These liners are intended for residual limbs with firm tissues and no bony prominences. We recommend a specific weight bearing socket with distal end contact in combination with a shuttle lock fitting (see OB shuttle locks 6A20=* and 6A30=*).

Important!

When using a locking liner the integrated matrix should end below the head of the fibula. Special orders with shorter matrix are possible.

4.2 Copolymer

These liners provide additional skin moisturising and can be used with a wide variety of limbs shapes. We recommend a total surface weight bearing socket in combination with an expulsion valve (see OB expulsion valves 4R140, 4R136, and 4R138=5).

Tip:

Thermal reshaping of copolymer liners. Pull the liner over a dry plaster positive. Wrap an elastic bandage around the unit if necessary.

6Y90/6Y92: Put this combination in an oven at a temperature of 70°C for two hours.

6Y93: Since it is made of a softer material, put this combination in an oven at 65°C for a maximum of 15 minutes.

4.3 Polyurethane

These liners are intended for sensitive residual limbs and/or with scars, grafts and/or bony prominences. We recommend a total surface weight bearing socket in combination with a Harmony® System (see 4R144, 4R150 or 4R148) or an expulsion valve (see OB expulsion valves 4R140, 4R136 and 4R138=5).

For additional information, please see our liner selection tool

5 Fabrication

We have created Technical Information guides (in multiple languages) to provide additional technical support for the different liners, as follows:

5.1 *Silicone, SiliconeGel and TechnoGel liners*

Technical Information 2.1.6

5.2 *TF Adapt*

Technical Information 2.3.1

5.3 *Copolymer and Polyurethane liners*

Technical Information 2.1.8 (this technical information is based for the Harmony® System. If using an expulsion valve – simply change the valve for the Harmony® component).

Please contact your local Otto Bock branch for your personal Technical Information copy or for additional information and support.

5.4 *Trimming a liner*

Attention!

During trimming, do not follow the edge of the socket but instead cut all around at an equal height slightly above / outside the edge of the socket.

Do not use wave edge scissors as this would result in a risk of the heavily prestressed liner tearing at the wave troughs of the cutting edge.

We recommend the following special scissors to shorten a liner: 719S7 and 717S20. These scissors are essential in order to cut the extremely durable textile material. With seamed textiles, use superglue to secure the proximal end of the seam. To trim a non-covered liner, we recommend using the trimmer (756L10) which will leave a bevelled edge.

The liner should not be deflected over the socket brim.

**Attention!**

If fitting with the Shuttle Lock, please make sure, that the patient has a socket with distal end contact. Furthermore the Shuttle Lock must always follow the long axis of the residual limb. Secure the pin screwed into the liner with Loctite® 242. Tighten it to **3 Nm** (2.2 lbf. in, 27 lbf. in).

6 Instruction for donning, doffing the liner and care

Please clean the liner daily by using running warm water and ph neutral soap. When using silicone liners, we recommend treating areas of the skin that will be subjected to a high amount of friction with 453H12 Derma Prevent before applying the liner. For more information please refer to the patient information.

**Attention!**

Powder of any kind (talcum, baby powder) between the residual limb and liner can cause skin irritations on the residual limb and signs of disintegration on the liner.

Please do not use any kind of powder on the inner side of the liner.

**Notice:**

Prior to first use, please wash the liner thoroughly.

7 Restrictions of Usage

**Attention!**

Use of a liner by another person may cause functional or hygienic risks.

The product is designed for use on **one patient only**. Re-use of the product is not permitted.

The daily duration of use and period of application are dependent on the medical indication.

8 Guarantee / Service Life

This product has a service life of 6 months defined by the manufacturer. During this time, material defects, manufacturing faults or premature wear are covered by the guarantee.

This guarantee applies, subject to normal wear and tear, when the products are used as intended, without unapproved modifications, following all Otto Bock instructions and requirements; and when they are fitted by or under the direct supervision of certified/ licensed practitioners who meet all Otto Bock product-specific training requirements, if any.

Otto Bock's sole obligation under this product guarantee shall be to repair or replace at no charge or to refund the cost of the product thereof, at Otto Bock's sole discretion.

Contact your local Otto Bock customer service representative for additional information.

9 CE Conformity

This product meets the requirements of the 93/42/EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH
U.S.A.–Minneapolis, Minnesota 55447